

# Bulk IgG's



# Bulk IgG'S

Immunoglobulins (IgG's) can be used as heterophilic blocking agents for many *in vitro* diagnostic applications to prevent false positive results. The addition of immunoglobulins are intended to block undesirable reactivity within the assay. A well-known example of heterophilic reactivity is HAMA. Mouse IgG has the capacity to reduce HAMA reactivity when testing human patient samples with mouse monoclonal antibodies.

## Specifications

**Purification:** Protein A or Protein G

**Purity:** >95% based on SDS PAGE

**Concentration:** 8-12 mg/ml

**Form:** Liquid or Lyophilized

**Unit Sizes:** 100 mg / 500 mg / 1 gram / 1 kg

Name	Part #
Rat IgG	Rt-003-Z
Mouse IgG	Mu-003-C
Human IgG	Hu-003-C
Goat IgG	Gt-003-Z
Rabbit IgG	Rb-003-C
Sheep IgG	Sh-003-Z
Chicken IgY	Ck-003-B