



# myReads® NGS Service Packages

Flexible, customizable options for different sample and project types

**Pick a package that fits your  
project needs or mix and match!**

Don't see anything you like?  
We offer à la carte services to meet your needs.

## STANDARD DNA

Appropriate for freshly-extracted, good quality DNA

- Includes QC, fragmentation, library preparation, 8-to-12-plex captures, and 1Gbp of PE150 AVITI sequencing per sample
- All samples must be high-molecular weight (>10kbp modal length), high quality (UV 260:280 1.7-1.9), and high mass (>500 ng of gDNA)
- Library preparation performance commitment\* for samples meeting our guidelines
- Complimentary QC report allowing you to make decisions about which samples to keep, drop, or supplement

## DEGRADED DNA

Appropriate for older specimens, herbarium, material, or samples of unknown quality

- Includes extended QC, fragmentation (where necessary), library preparation with specialized chemistry (where necessary), 4-to-6-plex captures, 1Gbp of PE150 AVITI sequencing per sample
- Samples may have a range of molecular weights, UV 260:280 values, and mass (>10 ng DNA)
- Performance is not guaranteed for degraded samples
- Complimentary QC report allowing you to make decisions about which samples to keep, drop, or supplement
- Perfect for samples for which you have no information or are unsure about gDNA quality/quantity!

## ANCIENT DNA - DOUBLE-STRANDED

Appropriate for museum, historical, or ancient DNA

- Includes QC, double-stranded library preparation with specialized chemistry (the "gold standard" for ancient DNA), 2-to-4-plex dual-round captures, 1Gbp of PE150 AVITI sequencing per sample
- All samples must be low-molecular weight (<1kbp modal length)
- Performance is not guaranteed for ancient samples
- Complimentary QC report allowing you to make decisions about which samples to keep, drop, or supplement

## ANCIENT DNA - SINGLE-STRANDED

Appropriate for museum, historical, or ancient DNA, cfDNA, single-stranded cDNA

- Includes QC, single-stranded library preparation with specialized chemistry, 2-to-4-plex dual-round captures, 1Gbp of PE150 AVITI sequencing per sample
- All samples must be low-molecular weight (<1kbp modal length)
- Performance is not guaranteed for ancient samples
- Complimentary QC report allowing you to make decisions about which samples to keep, drop, or supplement

## LONG-INSERT DNA

Appropriate for DNA or cDNA intended for PacBio sequencing

- Includes QC, long-insert precursor library preparation, 3-plex captures, 24 samples per SMRTbell preparation & PacBio Revio cell
- All samples must be high-molecular weight (>10kbp modal length), high quality (UV 260:280 1.7-1.9), and high mass (>5 µg DNA)
- Library preparation performance commitment\* for samples meeting our guidelines
- Complimentary QC report allowing you to make decisions about which samples to keep, drop, or supplement

## STANDARD RNA

Appropriate for high-quality RNA or mRNA

- Includes QC, DNAase treatment, purification, 8-to-12-plex captures, 1 Gbp of PE150 AVITI sequencing per sample
- All samples must be high quality (UV 260:280 1.9-2.1), and high mass ( $\geq$  250 ng RNA)
- Additional options for mRNA capture & rRNA depletion from human/mouse/rat samples
- Library preparation performance commitment\* for samples meeting our guidelines
- Complimentary QC report allowing you to make decisions about which samples to keep, drop, or supplement

### A LA CARTE SERVICES

*Your project, your way*

Need fewer samples per capture? DNA extraction?  
More sequencing? Bioinformatics? We've got you covered.

Talk to our dedicated team of scientists about the best custom options for your project today.

\* See *myReads Services Policies* document for details.  
Daicel and Daicel Arbor Biosciences are registered trademarks of Daicel Corporation.



web: [www.arborbiosci.com](http://www.arborbiosci.com)  
email: [info@arbor.daicel.com](mailto:info@arbor.daicel.com)  
phone: 1-734-998-0751  
twitter: @ArborBio



Download doc

**myReads**<sup>®</sup>  
Customized NGS Services  
[arborbiosci.com/myReads](http://arborbiosci.com/myReads)  
08/2023